

OXYGEN THERAPY

NONREBREATHING MASK

WITH SAFETY VENT

- high concentration
- half open, low-resistance check valve prevents rebreathing and allows exhaled gas to escape
- flow rate: 5-10 LPM
- clear, soft vinyl
- latex-free
- clean
- single use
- individually packed



ADULT NONREBREATHING ELONGATED MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
1059	with 210 cm of oxygen supply tubing	50
41069	without oxygen supply tubing	

PAEDIATRIC NONREBREATHING ELONGATED MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
41058	with 210 cm of oxygen supply tubing	50

NONREBREATHING MASK

- high concentration
- closed, low-resistance check valve prevents rebreathing and allows exhaled gas to escape
- flow rate: 10-12 LPM
- clear, soft vinyl
- latex-free
- clean
- single use
- individually packed



ADULT NONREBREATHING ELONGATED MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
41060	with 210 cm of oxygen supply tubing	50

MATERIALS

- mask: PVC with metal bracket
- tubing: PVC
- reservoir: PVC / 750 ml
- valve: silicone
- neck band: polychloroprene

SEE-THRU® MASK

- medium concentration
- flow rate: 5-10 LPM
- clear, soft vinyl
- latex-free
- clean
- single use
- individually packed



ADULT SEE-THRU MASK		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
41040	with 210 cm of oxygen supply tubing	50

PAEDIATRIC SEE-THRU MASK		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
41035	with 210 cm of oxygen supply tubing	50

SEE-THRU® MASK

- medium concentration
- flow rate: 5-10 LPM
- clear, soft vinyl
- latex-free
- clean
- single use
- individually packed



ADULT SEE-THRU ELONGATED MASK		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
1041	with 210 cm of oxygen supply tubing	50
1049	without oxygen supply tubing	

PAEDIATRIC SEE-THRU ELONGATED MASK		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
41042	with 210 cm of oxygen supply tubing	50
41050	without oxygen supply tubing	

SEE-THRU® MASK

- medium concentration
- flow rate: 5-10 LPM
- clear, soft vinyl
- adaptor swivels to accommodate patient position
- latex-free
- clean
- single use
- individually packed



ADULT SEE-THRU ELONGATED MASK		<i>HUDSON RCI</i>
REF.	DESCRIPTION	QTY
41007	with 210 cm of oxygen supply tubing & reservoir	50

MATERIALS

- mask: PVC with metal bracket
- tubing: PVC
- reservoir: PVC/750 ml
- neck band: polychloroprene

3-IN-1 MASK

3-in-1 masks can be used as a medium concentration, high concentration or nonrebreathing mask.

Each unit includes:
SEE-THRU mask with flapper valve,
reservoir and oxygen connector

- high concentration
- flow rate: 5-10 LPM
- clear, soft vinyl
- latex-free
- clean
- single use
- individually packed



ADULT 3-IN-1 ELONGATED MASK *HUDSON RCI*

REF.	DESCRIPTION	QTY
41061	with 210 cm oxygen supply tubing	50
41063	with 210 cm STAR LUMEN tubing (kink-resistant)	

MATERIALS

- mask: PVC with metal bracket
- tubing: PVC
- reservoir: PVC / 750 ml
- connector: polypropylene
- neck band: polychloroprene

MULTI-VENT® MASKS

VENTURI SYSTEMS

- features colour-coded, air-entrainment low and medium-concentration diluters
- includes adaptor for high humidity
- locking ring secures flow setting
- latex-free
- clean
- single use



Low-concentration diluter (green) for 24%, 26%, 28% and 30%



Medium-concentration diluter (white) for 35%, 40% and 50%

ADULT MULTI-VENT ELONGATED MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
41088	with 210 cm of oxygen supply tubing	50

PAEDIATRIC MULTI-VENT ELONGATED MASK

HUDSON RCI

REF.	DESCRIPTION	QTY
41089	with 210 cm of oxygen supply tubing	50

MULTI-VENT Total Gas Flow to Patient*

OXYGEN % SETTING	SUGG. OXYGEN FLOW	TOTAL GAS FLOW
24%	3 LPM	79 LPM
26%	3 LPM	47 LPM
28%	6 LPM	68 LPM
30%	6 LPM	53 LPM
35%	9 LPM	50 LPM
40%	12 LPM	50 LPM
50%	15 LPM	41 LPM


* At specific diluter and flow settings

SELECT-A-VENT® AIR ENTRAINMENT MASKS

VENTURI SYSTEMS

Available as a kit (including six diluters: 24%, 28%, 31%, 35%, 40% and 50%) or as a single diluter.

- complete with 210 cm of oxygen tubing supply
- includes adaptor for high humidity
- latex-free
- clean
- single use



MATERIALS

- mask: PVC with metal bracket
- tubing: PVC
- connector: polypropylene
- neck band: polychloroprene

ADULT SELECT-A-VENT AIR

HUDSON RCI

ENTRAINMENT ELONGATED MASKS

REF.	DESCRIPTION	QTY
41098	kit with 210 cm of oxygen supply tubing	50